OTFI-01X0 Fiber Illuminator

User's Manual ©

1. Scope

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3.2. WHEN CHANGING FIBER, DEVICE SHOULD BE TURNED OFF

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During operation, some internal components can become hot. To avoid the possibility of injury, always turn the device off and allow 30 seconds for cooling prior to changing the fiber or adapter or otherwise exposing the internal components.

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3.3. PROVIDE ADEQUATE VENTILATION AT ALL TIMES

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Allow clearance for the ventilation slots as described in section 1.1 to insure the device does not overheat and trigger the internal shut-off mechanism.

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3.4. DIRECT VIEWING OF EXITING LIGHT CAN BE HARMFUL

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Avoid direct viewing of the light. All LEDs and LED-based products are sources of light which can in some circumstances cause damage to the eye. The potential for any eye injury will depend upon many factors including but not limited to:

- duration of exposure and extent of pupil dilation prior to exposure
- distance from light source to eye
- drive current supplied to LEDs determined by the intensity light setting
- beam pattern (degree of focus)

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3.5. USE ONLY THE POWER CONVERTER SUPPLIED WITH THE PRODUCT

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Use of other power converters can damage product, will void warranty, and may result in unsafe operation

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3.6. PLUG AC CORD ONLY INTO SUITABLY GROUNDED AC RECEPTACLES

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3.7. ALWAYS USE THE DEVICE IN A MANNER SPECIFIED BY THE MANUFACTURER

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If the device is used in a manner not specified by the manufacture, the protection provided by the equipment may be impaired.

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3.8. DO NOT USE FOR MEDICAL APPLICATIONS

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This product is not intended for medical use.

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4. General Information

4.1. Intended use of equipment

The Excelitas Technologies OTFI-01X0 series Fiber Illuminator is a LED-based fiber optic illuminator designed for general microscopy and for industrial applications such as assembly illumination, machine vision, boroscopes, and inspection.

4.2. Technical specifications

Technical data	Value
Light output at maximum intensity, models OTFI-0100 & OTFI-0180	905 lumens typical ¹
Light output at maximum intensity, models OTFI-0110	710 lumens typical ¹
Color temperature	5,000 to 6,000°K ²
Input voltage	100 to 240V AC, 1.5 A, 50 to 60Hz
Overvoltage Category	II
Dimensions	3.9" W x 5.1" H x 7.0" L (9.9cm W x 13.0cm T x 17.8cm L)
Weight, models OTFI-0100 & OTFI-0110	5.9 lbs/2.7kg
Weight, model OTFI-0180	4.2lbs/1.9kg

Environmental use conditions	Value
Ambient temperature	0°C to 40°C
Relative humidity	80% max at 31°C decreasing linearly to 50% at 40°C
Altitude	Up to 2000 m/6600 ft
Pollution degree	2

¹Measured luminous flux directly out of Fiber Illuminator

²Measured CCT from the distal end of a 5 mm glass fiber bundle with NA = 0.6. CCT may vary if measured with optical fibers having different characteristics.

4.3. Instructions for use

See sections 5 (Equipment Installation) and sections 6 (Equipment Operation).

4.4. Electrical Interference

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the equipment and receiver
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected
- Consult the dealer or an experienced radio/TV technician for help

4.5. Equipment components

The product includes the following components:

- Fiber illuminator
- External power converter (AC power supply) with ferrite and attached cord
- AC power cord (sold separately)
- Adapter-holder for fiber optic bundles (sold separately)

5. Equipment installation

5.1. Positioning the device

Choose a suitable location that provides the following:

- Adequate clearance to connect the power converter per section 5.2.
- Within reach of a grounded AC outlet.
- Adequate clearance to operate the controls per sections 6.2 and 6.3.
- Clearance for ventilation per section 1.1.

5.2. Connections to the AC power supply

5.2.1. Insert the male plug from the external power converter into the female power jack on the back of the unit.


5.2.2. Plug the AC cord into the receptacle on the external power converter and into the AC outlet.

5.2.3. Plug AC cord only into suitably grounded AC receptacles

5.2.4. Use only the power converter provided with the product - use of other power converters can damage product, will void warranty, and may result in unsafe operation.

Inserting power converter plug into rear of Fiber Illuminator

AC plug connection

5.3. Ventilation requirements

Do not block the ventilation slots provided for the internal cooling fan. Slots are located on both sides of the unit and on the back.


6. Equipment operation

6.1. Attaching a fiber optic adapter

If fiber adapter is not already attached to Fiber Illuminator unit, orient the adapter such that the adapter code designation is viewable on the top of the adapter as shown below. Insert the adapter into the Fiber Illuminator aligning the three pins with the corresponding three slots, push and twist clockwise until the thumbscrew is oriented horizontally as shown.


6.2. Turning the device on and off

To turn the device on, push the power switch on the lower left of the front panel to the ON position. A green indicator light will turn on when the device is powered on.

To turn the device off, push the switch to the OFF position.

6.3. Changing intensity

To change the intensity of the light output, rotate the *Intensity* dial on the bottom right of the front panel. Rotating clockwise increases light output and rotating counterclockwise decreases intensity.

6.4. Instructions for cleaning

6.4.1. The exterior of the device can be gently cleaned with a soft cloth dampened with water and mild detergent. Do not use any solvents or cleaner containing chlorine as these will damage the device and may void your warranty.

6.4.2. The ventilation slots can be cleaned by gently blowing aerosol- based dusting air into the slots. Do not blow compressed "shop-air" into the device as such air often contains oil and water mist which may degrade the internal optics.

7. Applicable standards

Category	Regulatory body	Title of standard	Standard
EMI	European Union	Electromagnetic Compatibility	Directive 2004/108/EC
Safety	European Union	Low Voltage Directive (LVD)	Directive 2006/95/EC
RoHS	European Union	Restriction of the use of Hazardous Substances in Electrical and Electronic Equipment	Directive 2002/95/EC
Immunity	EN	Electrical equipment for measurement, control and laboratory use- EMC requirements- Part 1: General requirements	EN 61326-1:2006
Emissions	EN	Industrial, Scientific and medical (ISM) Radio Frequency Equipment- Electromagnetic Disturbance Characteristics- Limits and Methods of Measurement	EN 55011:2007 Amendment A2:2007, Group 1, Class B
Emissions	EN	Electromagnetic compatibility (EMC)-Part 3-2: Limits- Limits for harmonic current emission (equipment input current ≤ 16 A per phase)	EN 61000-3-2:2006
Emissions	EN	Electromagnetic compatibility (EMC) - Part 3-3: Limits - Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current ≤ 16 A per phase and not subject to conditional connection.	EN 61000-3-3:1995 Amendment A1:2001 to EN 61000-3-3:1995 and Amendment A2:2005 to EN 61000-3-3:1995
Emissions	ICES	Spectrum Management and Telecommunications Policy Interference-Causing Equipment Standard Digital Apparatus	ICES-003 Issue 4 Class B February 2004
Emissions	FCC	FCC Rules and Regulations/Unintentional Radiators Class B Digital Devices	FCC Part 15, Subpart B, Sections 15.107a and 15.109a Class B Digital Devices
Safety	IEC/EN	Safety requirements for electrical equipment for measurement, control, and laboratory use. Part 1- General Requirements	IEC/EN 61010-1:2001
Safety	UL	Safety requirements for electrical equipment for measurement, control, and laboratory use. Part 1- General Requirements	UL 61010-1: 2005-07-22
Safety	CSA	Safety requirements for electrical equipment for measurement, control, and laboratory use. Part 1- General Requirements	CSA C22.2 NO 61010-1-2004-07

8. Markings

9. Patents

The OTFI-01XX products are covered by "U.S. Pat. No. D629,947". Other patents pending.

 Wheeling, IL 60090 USA www.excelitas.com	Title OTFI-01X0 Fiber Illuminator User's Manual		Number OTFI-01X0-17	
	Author J. Swayne	Date 1/13/2012	Revision C	Page 9 of 10

10. Limited Warranty

Excelitas Technologies warrants the OTFI-01XX product to be free from defects in material and workmanship and to be in conformance with the written specification for a period of 12 months from date of purchase.

If any defect in material or workmanship or failure to conform to such specification is found, the Purchaser should promptly notify Excelitas Technologies. After a Returned Material Authorization number is assigned by Excelitas Technologies, Purchaser may return the product to Excelitas Technologies, carrying charges prepaid. Products will not be accepted for repair, replacement, credit or refund, without the written authorization of and in accordance with Excelitas Technologies instructions.

Excelitas Technologies shall analyze the failures, making use, when appropriate of technical information provided by Purchaser relating to the circumstances surrounding the failures. At Excelitas Technologies option, we will repair or replace the product found to be defective, and shall return the product carrying charges prepaid. Excelitas Technologies correction of any defects by the grant of credit, replacement or repair shall constitute fulfillment of all of its obligations and liability to the Purchaser hereunder.

Excelitas Technologies is not responsible for damage to its product caused by misuse, neglect, accident, shipping, abuse, maintenance, cleaning procedures, use or attempts to operate above its rated capacity intentionally or that otherwise deviate from the parameters established in the user integration manual and applicable specifications; or to products that have been improperly installed, stored, maintained, repaired or altered by anyone other than Excelitas Technologies; or have had their serial numbers or month and year of manufacture or shipment removed, defected or altered. Any actions cited above shall terminate this Warranty and shall relieve Excelitas Technologies from any further responsibility.

Excelitas Technologies shall not be liable for any incidental, special, or consequential damages in any claim action, suit or proceeding arising under this Warranty or any other part of the agreement of sale between Excelitas Technologies and the Purchaser, nor shall there be any liability hereunder for claims for labor, loss of profits or good will, repairs or other expenses incidental to replacement.

This warranty does not extend to any system into which a product is incorporated. No other warranty, including warranties of merchantability or fitness for a particular purpose is given with respect to such service or any other service provided by Excelitas Technologies under this Agreement. This warranty applies only to Purchaser and may not be assigned or extended by Purchaser to any of its customers or other users of the Items. Excelitas Technologies will not accept any returns from Purchaser's customers or users of Purchaser's products.

Product specifications are subject to change without notice.

Service Contact Information

For service please contact:

Excelitas Technologies LED Solutions
160 E. Marquardt Drive
Wheeling, IL 60090
USA
T: (847) 537- 4277
F: (847) 537- 4785
E-mail: ledsolutions.na@excelitas.com
www.Excelitas.com

© Copyright USA

 Wheeling, IL 60090 USA www.excelitas.com	Title OTFI-01X0 Fiber Illuminator User's Manual		Number OTFI-01X0-17	
	Author J. Swayne	Date 1/13/2012	Revision C	Page 10 of 10